




Crystals

for X-Ray
spectrometry



The spectral analysis of X-rays emitted by a sample after irradiation is both a powerful qualitative and quantitative analytical technique. It is based on the following phenomenon: an atom relaxes after excitation by emitting X-ray radiations at specific wavelengths, which reveal the identity of the emitting species.


For this spectral analysis Saint-Gobain Crystals supplies two key components:

- The monochromating crystals
- The scintillation detectors

An X-ray spectrometer basically consists of:

- An excitation source which may be
 - Either a primary X radiation, in which case one refers to X-ray fluorescence spectrometry.
 - Or an electron beam, inducing a so-called direct emission, used in microprobes and scanning electron microscopes.
- A monochromating crystal which is used to disperse the various spectral components of the incident beam.
- A detector in order to measure the intensity of the various spectral lines as singled out by the monochromator.

The detector offered by Saint-Gobain Crystals is a X Scintibloc® which combines a NaI(Tl) scintillator directly coupled to a photomultiplier with a low absorbing MIB or beryllium entrance window (see brochure "Scintillation detectors").



Monochromating crystals

A monochromating crystal behaves in X-ray spectrometry as does a diffraction grating in optics. When rotated with respect to the incident polychromatic beam (see figure), it will diffract the spectral component along with direction to satisfy Bragg's law, namely:

$$2d \sin \theta = n \lambda$$

where integer n refers to the diffraction order.

Hence, the most important characteristic of a monochromating crystal is the double atomic spacing $2d$, which gives the largest wavelength to be diffracted.

The range of monochromators supplied by Saint-Gobain Crystals can be found in the table (next page), along with the usual surface finish, within our control means, to the best intensity-resolution compromise. The optimum depends on each specific case and strongly reflects the nature of the set-up.

Products available

Monochromating crystals can be supplied in the two following shapes:

- Flat
- Curved onto a holder

Flat plates can be supplied unmounted or mounted into holders suitable for industrial X-ray fluorescence spectrometers. Other types of holders may also be supplied on request. The standard orientation accuracy provided is 10 minutes. On special request a one-minute accuracy can be ensured.


The curved plates used in such instruments as microprobes and scanning electron microscopes are always supplied on tailor-made holders.

Two main types of focusing configurations may be considered:

- The Johann geometry
- The Johansson geometry

3D curved crystal optics can be provided according to your special finishing requirements:


- spherical
- toridal
- ellipsoidal
- conical
- others (please inquire)



TIAP crystal



Natural quartz



PET crystal



Natural beryl


Catalogue of Saint-Gobain Crystals monochromators

Crystal	Lithium fluoride (1)			Quartz (1) (2)		Indium Antimonide	Silicon (1) (2)		Germanium (1) (2)	
	Chemical formula	LiF			SiO ₂		InSb	Si		Ge
Crystal system	Cubic			Hexogonal (3)		Cubic	Cubic		Cubic	
Parameters:	4.027			4.913 4.913 5.405		6.48	5.431		5.658 8.74	
Reflecting planes orientations	(200)	(220)	(420)	(10 $\bar{1}1$)	(10 $\bar{1}0$)	(111)	(111)	(220)	(111)	(220)
2d in Å	4.027	2.848	1.801	6.684	8.514	7.480	6.271	3.840	6.532	4.000
Usual surface finish	Cleaved or	Treated	Treated	Polished	Polished	Polished	Polished	Polished	Polished	Polished
Reflectivity	Intense	Intense	Average	Good	Good	Intense	Intense	Average	Intense	Intense
Calibration elements	Mo, Fe, Ti	Mo, Fe	Mo	Cu	Cu	Si	Cu	Cu	Cu	Cu
Common applications	From K to heavy elements Heavy elements Lines splitting			As Ge (111) As PET		Quantitative analysis of silicon	Extinction of even order spectral lines		Extinction of even order spectral lines	

Crystal	Pentaerythritol PET	Ethylene Diamine Dextratrate EDDT	Ammonium Dihydrogen Phosphate ADP	Beryl (2)	Acid Phthalates		
					Thallium TIAP	Rubidium RbAP	Potassium KAP
Chemical formula	C(CH ₂ OH) ₄	C ₆ H ₁₄ N ₂ O ₆	NH ₄ H ₂ PO ₄	3BeO·Al ₂ O ₃ ·6SiO ₂	CO ₂ HC ₆ H ₄ CO ₂ TI	CO ₂ HC ₆ H ₄ CO ₂ Rb	CO ₂ HC ₆ H ₄ CO ₂ K
Crystal system	Quadratic	Monoclinic	Quadratic	Hexagonal	Orthorhombic	Orthorhombic	Orthorhombic
Parameters:	6.16 8.74 5.959 105° 33'	8.97 8.808 7.542	7.530 7.530 9.17	9.21 9.21 12.95	6.63 10.54 13.06	6.55 10.02 13.32	6.46 9.61
Reflecting planes orientations	(002)	(020)	(101)	(10 $\bar{1}0$)	(001)	(001)	(001)
2d in Å	8.740	8.808	10.648	15.950	25.900	26.120	26.640
Usual surface finish	Cleaved or treated	Polished	Polished or treated	Polished	Cleaved	Cleaved	Cleaved
Reflectivity	Intense	Average	Average	Average	Intense	Intense	Good
Calibration elements	Al, Si	Al	Mg	Mg	Na, Mg	Na	Na
Common applications			Mg	Na and following elements	F to Al	Na to Al, up to F in emission probes	Na to Al, up to F in emission probes

- (1) Other cuts can be supplied upon request.
- (2) Noteworthy for its high radiation damage threshold (synchrotron).
- (3) The quartz unit cell is rhomboedric. The parameters given here correspond to a multiple hexagonal cell, defined as for beryl by $a = b = c$
 $\alpha = \beta = 90^\circ$ $\gamma = 120^\circ$

May also be supplied on request:
 Topaz (303) 2d = 2.712 Å
 Mica (002) 2d = 19.800 Å
 CsAP (001) 2d = 25.650 Å
 THM (010) 2d = 52.500 Å



Saint-Gobain Crystals will choose the most appropriate technique according to the type of crystal, its dimensions and the radius of Rowland circle to be achieved.

The orientation accuracy is better than 10 minutes, except for the Johann cleaved configuration, where it is 1 minute.


Other types of curvature may be investigated on request. For example, curvatures on holders shaped as logarithmic spiral, elliptic, parabolic, even spheric designs, interesting in plasma or synchrotron radiation study and astrophysics.

Manufacturing capabilities will strongly depend upon the crystal nature and dimensions as well as on the curvature radii.

The Johann geometry

A thin plate, produced by one of the two following methods:

- Cleavage for LiF (200)-PET-TIAP-RbAP-KAP
- Machining for other materials, is further cylindrically curved and glued upon a holder of curvature radius $2R$. It is possible to show that a beam emitted by a source at S is approximately focused at F . The source and focus are both located on the so-called Rowland circle whose radius is R




The Johansson geometry

Two different types of Johansson configurations, theoretically leading to perfect focusing, are considered:

- Single machining Johansson
- Double machining Johansson

The illustrations represent the fabrication steps of Johansson plates in both techniques.



B.P. 521
77794 Nemours Cedex, France

Telephone:
USA: 1-800-830-3594
Asia: +81(0)3 3263-0559
Europe: +33 1(0) 64 45 10 15

www.photonic.saint-gobain.com



Natural beryl

The data presented is believed to be correct but not guaranteed to be so.

©2004 Saint-Gobain Ceramics